



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Micro Centrifuge Bench-Top

**Product Code :**  
LBNY-0001-830008



**Description :**

Micro Centrifuge Bench-Top

**Technical Specification :**

These centrifuges can be recognised by the red cover of the centrifuges as well the rotors which also have red covers.

The table-top and floor standing centrifuges have an identical drive and mostly all of the rotors of the can fit in every centrifuge of the range.

The models differ in the sizes of the chamber and the power of the machine.

Therefore the maximum speed differs across the range.

Cooling: non refrigerated

Maximum capacity: 24x1.5/2.0ml

Speed: 100-15060 rpm (RCF: 21300 xg)

Maximum RCF with fixed angle rotor 21300 xg

Acceleration time: 15 seconds

Timer: 10 seconds to 9.59 hours, with continuous run function

Noise level: 75 dB(A)

---

# Laboratory Instrument India

