




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Fluctuating Pictures Projection Demonstrator

**Product Code :**  
LBNY-0004-101000417



**Description :**

It is used to demonstrate dynamic automatic cycle display of transverse wave propagation; dynamic display cycles of transverse wave propagation; dynamic display of particle movement trend away from its balance place; static display of every frame waveform image and static double wave display.

**Technical Specification :**

Dimension (L\*W\*H CM): 53.3\*23.4\*6



**Laboratory Instrument India**