



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bench-Top Refrigerated Centrifuge

Product Code :
LBNY-0001-830009



Description :

Bench-Top Refrigerated Centrifuge

Technical Specification :

Brushless DC frequency motor, longer life and low noise.

Flexible axle driven system which drive the rotor directly, smooth in operation.

Imported compressors, fluorine-free, double cycle cooling, cold and hot alternating easily, environment pollution-free and precise in temperature control.

Microprocessor control speed, temperature and time, digital display, 10 kinds of accelerating and decelerating speed.

Temperature range -9 to 40oC

Maximum speed 4000 rpm

Digital display

Timer, programmable up to 99min or continuous

Fixed-angle rotor, 6x250ml

Swing-out rotor, 4x100ml Accessories

Rotor adaptor; 8x15ml tubes, 6x15ml tubes

Laboratory Instrument India

