



Email : sales@laboratoryinstrumentindia.com

Product Name :
Floor Refrigerated Centrifuge

Product Code :
LBNY-0001-830010



Description :

Floor Refrigerated Centrifuge

Technical Specification :

Microcomputer controlled programmable operation with touch panel
TFT true-colour LCD wide screen displays set and run conditions for easy monitoring
Brushless DC Motor ensures low maintenance, high efficiency, constant speed of rotation
Fully stainless steel structure and a stainless steel cavity makes the machine safe and more efficient
Temperature range -9 to 40oC
Maximum speed 4000 rpm
Digital display
Timer, programmable up to 99min or continuous
Fixed-angle rotor
Swing-out rotor
Accessories



Laboratory Instrument India

