



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Vortex and/or Multivortex

**Product Code :**  
LBNY-0001-830011



**Description :**

Vortex and/or Multivortex

**Technical Specification :**

Intended for intensive stirring of bacterial and yeast cell, washing from the culture medium and extraction of metabolites and enzymes from cells and cell cultures.

It is different from by the possibility of mixing up to 32 tubes simultaneously.

32-socket universal platform for Eppendorf type tubes up to 1.5 ml (1.5/0.5/0.2 ml - 16/8/8 sockets)

Platform for mixing single tube up to 15 ml.

6-socket platform PV-6/10 for 10 ml tubes (max. tube diameter 15 mm)

Model of shaker: orbital shaker vortex

Shaking weight: 70 g

Shaker speed range: 3000 rpm



**Laboratory Instrument India**