



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Wind Cup Type Wind Speed Meter

**Product Code :**  
LBNY-0006-10100032



**Description :**

It can measure the instantaneous wind speed, average wind speed, instantaneous wind level, average wind speed level, the corresponding wave higher parameters and data latch function.  
Composed of the wind speed sensor and the hostpanel with LCD display and function keys.

**Technical Specification :**

Wind Cup Type Wind Speed Meter



**Laboratory Instrument India**