



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Minimum Thermometer

**Product Code :**  
LBNY-0006-10100034



**Description :**

Measuring range:  $-52^{\circ}\text{C}$  to  $+41^{\circ}\text{C}$ , accuracy:  $\pm 1^{\circ}\text{C}$ , resolution:  $1^{\circ}\text{C}$ .

**Technical Specification :**

Minimum Thermometer



**Laboratory Instrument India**