




Email : sales@laboratoryinstrumentindia.com

Product Name :
6 Hole Test Tube Rack

Product Code :
LBNY-0006-10100039



Description :

The six hole type test tube stand is made up of plastic and the diameter is 20 ± 1 mm, the distance of each hole center is 25 ± 1 mm.

The distance between the orifice plate and the base is 85 ± 2 mm.

It has six tube columns, of which the diameter is 10 ± 1 mm and the length is about 65 ± 2 mm.

Technical Specification :

6 Hole Test Tube Rack



Laboratory Instrument India