



Email : sales@laboratoryinstrumentindia.com

Product Name :
Towel UV Sterilizer

Product Code :
LBNY-0001-830015



Description :

Towel UV Sterilizer

Technical Specification :

Built-in UV lamp to sterilize different types of tools.

Three timer setting modes (30-60-90min) for your selection LED time display panel to precisely calculate the working time.

Glass magnetic door for prevention of the UV ray and viewing of sterilizing process.

Auto-shutoff design: UV lamp will go off automatically when the draw is open.

8L capacity

Built-in UV lamp to sterilize different types of tools

Multiple timer setting modes

LED time display panel to precisely calculate the working time

Glass magnetic door for prevention of the UV ray and viewing of sterilizing process

Auto-shutoff design: UV lamp will go off automatically when the draw is open

Fuse protected, safety to use



Laboratory Instrument India

