



Email : sales@laboratoryinstrumentindia.com

Product Name :
Resistance Decade Box

Product Code :
LBNY-0002-620013



Description :

Resistance Decade Box

Technical Specification :

For the completion of our sales program in decade with shielding plastic housing for education, maintenance, trouble-shooting, working standard research, design and developing.
The rotary switches enables precise, decadic adjustment of the resistance-, capacitance- or inductance values.
Variable from 1 to 10 m² Resistance

Resistance Decade Box, Resistance Decade Box Equipments, Resistance Decade Box Tools, Resistance Decade Box Tool Kits, Resistance Decade Box Manufacturers, Resistance Decade Box Suppliers from India, China, Kenya



Laboratory Instrument India