



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Resistance Decade Box

**Product Code :**  
LBNY-0002-620013



**Description :**

Resistance Decade Box

**Technical Specification :**

For the completion of our sales program in decade with shielding plastic housing for education, maintenance, trouble-shooting, working standard research, design and developing.  
The rotary switches enables precise, decadic adjustment of the resistance-, capacitance- or inductance values.  
Variable from 1 to 10 m<sup>9</sup> Resistance

Resistance Decade Box, Resistance Decade Box Equipments, Resistance Decade Box Tools, Resistance Decade Box Tool Kits, Resistance Decade Box Manufacturers, Resistance Decade Box Suppliers from India, China, Kenya



**Laboratory Instrument India**