



Email : sales@laboratoryinstrumentindia.com

Product Name :
UV Cabinet PCR Workstation

Product Code :
LBNY-0001-830016



Description :

UV Cabinet PCR Workstation

Technical Specification :

2 built-in shortwave (254 nm) UV tubes for decontamination between experiments
Timer sets UV exposure up to 12 h
Safety shut-off switch automatically turns the UV light off when door is opened
Keylock prevents accidental exposure of samples to UV
Unique, easy-clean antimicrobial coating on the stainless steel and aluminum surfaces
Hinged door flips up for easy access to the work area
Built-in power outlets for operation of equipment inside the work area
254nm UV radiation
Anti-microbial coating on stainless steel and aluminium surfaces
Automatic safety shutoff switch
Blocks all UV light below 400nm
Touch screen interface to control all system functions
HEPA filter system

Laboratory Instrument India

