



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
UV Cabinet PCR Workstation

**Product Code :**  
LBNY-0001-830016



**Description :**

UV Cabinet PCR Workstation

**Technical Specification :**

2 built-in shortwave (254 nm) UV tubes for decontamination between experiments  
Timer sets UV exposure up to 12 h  
Safety shut-off switch automatically turns the UV light off when door is opened  
Keylock prevents accidental exposure of samples to UV  
Unique, easy-clean antimicrobial coating on the stainless steel and aluminum surfaces  
Hinged door flips up for easy access to the work area  
Built-in power outlets for operation of equipment inside the work area  
254nm UV radiation  
Anti-microbial coating on stainless steel and aluminium surfaces  
Automatic safety shutoff switch  
Blocks all UV light below 400nm  
Touch screen interface to control all system functions  
HEPA filter system

---

# Laboratory Instrument India

