



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Measuring Line

**Product Code :**  
LBNY-0006-101000109



**Description :**

The dividing mark is clear, should be combined with the test rope. Both ends were combined tightly with plastic ring.

Polyester fiber wire, measuring range: 50m, thickness:  $\approx 2.5\text{mm}$ , resolution: 1m, accuracy:  $\pm 10\text{mm}$ . Total accuracy:  $\pm 250\text{mm}$ .

**Technical Specification :**

Measuring Line



**Laboratory Instrument India**