



Email : sales@laboratoryinstrumentindia.com

Product Name :
Incubators for Microbiology

Product Code :
LBNY-0001-830021



Description :

Incubators for Microbiology

Technical Specification :

Triple walled construction, the incubator has its inner chamber made from stainless steel sheet, center with G.I. sheet and outer cabinet is from mild steel sheet.

Double door inner is full view thick plate glass and outer is doubled walled insulated.

Supply is complete with plug, cable, neon, switches etc and G.I. wire netting shelves.

Application: Bacterial and microbiology research

Operating temperature range 5 to 75°C

Spatial temperature deviation at 37°C; $\pm 0.6^\circ\text{C}$

Chamber volume 150-200L

Number of shelves; 15-20

Maximum shelf load (kg); 25-50



Laboratory Instrument India

