



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
The Circular Area of PI Demonstrator

**Product Code :**  
LBNY-0006-101000126



**Description :**

The circumference ratio demonstration consists of a circle wheel ( 100mm diameter ) and a scale below the tray. The circle area demonstration consists of a 200mm diameter circle divided into 16 equal parts of the fan-shaped small pieces, and each one be divided into two small pieces .

**Technical Specification :**

Dimension (L\*W\*H CM): 46\*36\*1.1



**Laboratory Instrument India**