



Email : sales@laboratoryinstrumentindia.com

Product Name :
Water Baths

Product Code :
LBNY-0001-830022



Description :

Water Baths

Technical Specification :

Displays actual and set point temperature simultaneously
Programmable pre-sets for frequently used temperatures
Integrated timer
Hinged gable cover
Calibration off-set feature
Reservoir made of 304 Stainless Steel
Reservoir drain (10 liter and larger baths)
Primary and automatic safety thermostats
Capacity; 20-30 L
Digital Display - Full-color TFT
Working Temperature Range °C - Ambient +5° to 100°
Temperature Stability °C - $\pm 0.1^\circ$
Temperature Calibration Capability - 1-point
Integrated timer
Hinged gable cover
Tank Material; Stainless Steel



Laboratory Instrument India