



Email : sales@laboratoryinstrumentindia.com

Product Name :
Laboratory Drying Oven Forced Convection

Product Code :
LBNY-0001-830023



Description :

Laboratory Drying Oven Forced Convection

Technical Specification :

The interior is manufactured from stainless steel and the exterior is robustly constructed from sheet steel finished in powder coated paint.

Laboratory Ovens are fitted with direct reading digital control unit and equipped with an over temperature protection analogue thermostat.

120 liters laboratory ovens are supplied complete with 2 shelves.

Internal volume; 120 L

Max. temperature; 300°C

Temperature stability at 150°C; $\pm 0,3^{\circ}\text{C}$

Temperature homogeneity; $\pm 2,0\%$

Timer; 0 to 9 999 min. (possible continuous operation)

Display; backlit colour LCD

Display resolution; $0,1^{\circ}\text{C}$

Number of shelves; 7

Max. shelf load; 20kg

Laboratory Instrument India

