



Email : sales@laboratoryinstrumentindia.com

Product Name :
Lab Orbital Shaker

Product Code :
LBNY-0001-830024



Description :

Lab Orbital Shaker

Technical Specification :

Three eccentric shaft balancing drive ensures shaking with uniform speed.
Speed adjustment settings permit both gentle and vigorous shaking.
Varieties of platforms are available for different glassware and vessels.
Speed: 10 – 300rpm including stir reverse
Platform Dimensions: (D x W) 400-600 x 400-600mm
Timer: 1min-99 hours w/Alarm or continuous
Orbit 35mm
Maximum Load 12 Kg
Selectable Ramp Rates 10 - 100rpm/sec



Laboratory Instrument India