



Email : sales@laboratoryinstrumentindia.com

Product Name :
Humidity Chamber

Product Code :
LBNY-0001-830025



Description :

Humidity Chamber

Technical Specification :

The humidity is controlled by humidistat from atmospheric humidity to 95% with a sensitivity of $\pm 3\%$.
Double walled inside anodized aluminium/stainless steel and outside mild steel painted in epoxy powder coating.
The door has glass Window. To work on 220/230 volts.A.C. inside water tank is made of stainless steel.
Relative humidity and temperature control chamber adjustable; 20 to 150° and 20 to 98% relative humidity
Indoor vertical stand
Internal Dimension; W*H*D (mm); 400*500*400
Power supply 220V/50Hz



Laboratory Instrument India