



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Humidity Chamber

**Product Code :**  
LBNY-0001-830025



**Description :**

Humidity Chamber

**Technical Specification :**

The humidity is controlled by humidistat from atmospheric humidity to 95% with a sensitivity of  $\pm 3\%$ .  
Double walled inside anodized aluminium/stainless steel and outside mild steel painted in epoxy powder coating.  
The door has glass Window. To work on 220/230 volts.A.C. inside water tank is made of stainless steel.  
Relative humidity and temperature control chamber adjustable; 20 to 150° and 20 to 98% relative humidity  
Indoor vertical stand  
Internal Dimension; W\*H\*D (mm); 400\*500\*400  
Power supply 220V/50Hz



**Laboratory Instrument India**