



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
BOD Meter

**Product Code :**  
LBNY-0001-830028



**Description :**

BOD Meter

**Technical Specification :**

A state of art, microprocessor based, dissolved oxygen instrument with many automated and specific application features, designed to provide laboratory results and accuracy under harsh industrial conditions.

All measurements are automatically compensated for temperature.

Manual salinity compensation in water allows direct determination of dissolved oxygen in saline waters.

With its internal barometer, the instrument is able to automatically compensate for changes in barometric pressure so there is no need for charts, altitude information or external barometric pressure information.

Resolution; 0.01 ppm

Saturation; 0.1% 0 - 40 mg/L DO range

Power: 12V, AC to DC adapter

DO probe.



**Laboratory Instrument India**

