



Email : sales@laboratoryinstrumentindia.com

Product Name :
Alcohol Burner

Product Code :
LBNY-0007-10100033



Description :

Alcohol burner, 150ml. Full set: with plastic lamp cap, ceramic core seat, cotton rope wick.
Style: the main body is composed of a thick glass bottle bottom small shoulder, for storage of fuel alcohol, the shoulder is not easy to be overturned, and can ensure the accidentally knocked over tilt when alcohol is not easy to overflow.

Technical Specification :

Alcohol Burner



Laboratory Instrument India