



Email : sales@laboratoryinstrumentindia.com

Product Name :
Data Logger

Product Code :
LBNY-0001-830030



Description :

Data Logger

Technical Specification :

Plug and Read: No cable or reader required, the USB connector enables quick connection to PCs to download data.

PDF Report: Auto encrypted and generated PDF report with password protection function; Personalized settings and management provided by free Elitechlog software.

Recorder measures and stores up to 8000 sets of high- resolution humidity

Temperature measurement ranges from -40°C to +85°C

Relative humidity; 0-100%..



Laboratory Instrument India