



Email : sales@laboratoryinstrumentindia.com

Product Name :
The Solution Conductivity Demonstration Device

Product Code :
LBNY-0007-10100042



Description :

The product is mainly composed of a demonstration board 1 sets, 5 sets of electrolyzer, lamp etc. Products use the power supply voltage: 4V ~ 6V. The demonstration board seat is made of ABS engineering plastic, which is provided with a socket 5 and socket, wiring terminal, wiring diagram printed on the demonstration panel, a demonstration board seat size 360 x 140 x 260mm; electrolytic cell by a transparent box, a box cover and a connecting wire electrode, etc.

Technical Specification :

The Solution Conductivity Demonstration Device



Laboratory Instrument India