



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Portal Ethylene Analyser

**Product Code :**  
LBNY-0001-830031



**Description :**

Portal Ethylene Analyser

**Technical Specification :**

Monitor controlled atmosphere and storage environments to ensure ethylene levels do not exceed critical thresholds.

Flux and headspace analysis: Collect accurate ethylene emission levels across treatments of single fruits or small-volume samples with low ethylene emission rates.

Digital display

Measuring from 0-200ppm

Lower detection limit of 0.025ppm.



**Laboratory Instrument India**