



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sedimentator

Product Code :
LBNY-0007-10100049



Description :

The centrifugal tube is made of glass, the volume is not less than 10ml, the centrifugal tube shape dimension is about $\Phi 21 * 105\text{mm}$; the centrifugal pipe frame is made of cold rolled sheet with thickness not smaller than 1.8mm, and the support is not less than 2 centrifugal tubes. The utility model relates to a hand-operated centrifugal sedimentation device, which is mainly composed of a transmission device, a centrifugal pipe, a centrifugal pipe sleeve, a centrifugal pipe rack and a table clamp.

Technical Specification :

Sedimentator



Laboratory Instrument India