



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Sedimentator

**Product Code :**  
LBNY-0007-10100049



**Description :**

The centrifugal tube is made of glass, the volume is not less than 10ml, the centrifugal tube shape dimension is about  $\Phi 21 * 105\text{mm}$ ; the centrifugal pipe frame is made of cold rolled sheet with thickness not smaller than 1.8mm, and the support is not less than 2 centrifugal tubes. The utility model relates to a hand-operated centrifugal sedimentation device, which is mainly composed of a transmission device, a centrifugal pipe, a centrifugal pipe sleeve, a centrifugal pipe rack and a table clamp.

**Technical Specification :**

Sedimentator



**Laboratory Instrument India**