



Email : sales@laboratoryinstrumentindia.com

Product Name :
Electric Drill

Product Code :
LBNY-0007-10100050



Description :

Main technical specifications: the maximum clamping diameter is 10mm; the no-load speed is 0~600 RPM / min. the battery voltage is 7.2V, 1.5Ah, charger charging voltage is 7.2V to 12V. The overall size of the suitcase is 350 x 260 * 90mm. The product is a rechargeable electric drill, which consists of an electric drill, a rechargeable battery (2 pieces) and a charger.

Technical Specification :

Electric Drill



Laboratory Instrument India