



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Water Electrolysis Experiment Device

**Product Code :**  
LBNY-0007-10100051



**Description :**

The seat is provided with a connecting device and a wiring polarity identification; electrode using carbon rod production, diameter 4mm, length 45mm, the closed room size was 120 mm \* 50 mm \* 25MM, volume 15ml; discharge nozzle is made of glass, the length of about 50mm, and the connecting hose products using milk hose.

The product is made of transparent polycarbonate material. It is made up of experiment base, electrode, airtight chamber, observation tube, air vent and funnel.

The overall size of the instrument is about 150 mm \* 60 mm \* 190mm, and the operating voltage of the instrument is DC 12V.

**Technical Specification :**

Water Electrolysis Experiment Device



**Laboratory Instrument India**