



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Lab Stereo Microscope

**Product Code :**  
LBNY-0001-830034



**Description :**

Lab Stereo Microscope

**Technical Specification :**

Fully coated optical system which provides a sharp and clear imaging with extreme flatness and contrast. This stereo microscope offers ergonomic features such as synchronized bi-lateral inter-pupillary distance adjustments.

Designed to be outfitted with a variety of lens accessories to increase or reduce magnification.

With a large zoom range and long-working distance, it can produce the best reproduction and true color imaging. Ideal for high-magnification, high-resolution observation in biological research, semi-conductor and other precision-oriented industries.

Head: Binocular (45° Inclined, 360° Rotatable)

Eyepiece: WF 10X, 20mm

Zoom Range: 0.75X-7.5X

Range of Magnification: 7.5X-75X

Working Distance: 102mm

Visual Field: 26.7mm - 2.7mm

Interpupillary Distance: 50mm - 75mm, Diopters regulation 5

Stand: Guide-Rail Stand

Illumination: Top & Bottom LED illumination.



**Laboratory Instrument India**