




Email : sales@laboratoryinstrumentindia.com

Product Name :
Inductance Decade Box

Product Code :
LBNY-0002-620015



Description :

Inductance Decade Box

Technical Specification :

For the completion of our sales program in decade with shielding plastic housing for education, maintenance, trouble-shooting, working standard research, design and developing.
The rotary switches enables precise, decadic adjustment of the resistance-, capacitance- or inductance values.
0.01mH to 10H Inductance

Inductance Decade Box, Inductance Decade Box Equipments, Inductance Decade Box Tools, Inductance Decade Box Tool Kits, Inductance Decade Box Manufacturers, Inductance Decade Box Suppliers from India, China, Kenya



Laboratory Instrument India