




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Inductance Decade Box

**Product Code :**  
LBNY-0002-620015



**Description :**

Inductance Decade Box

**Technical Specification :**

For the completion of our sales program in decade with shielding plastic housing for education, maintenance, trouble-shooting, working standard research, design and developing.  
The rotary switches enables precise, decadic adjustment of the resistance-, capacitance- or inductance values. 0.01mH to 10H Inductance

Inductance Decade Box, Inductance Decade Box Equipments, Inductance Decade Box Tools, Inductance Decade Box Tool Kits, Inductance Decade Box Manufacturers, Inductance Decade Box Suppliers from India, China, Kenya



**Laboratory Instrument India**