



Email : sales@laboratoryinstrumentindia.com

Product Name :
COD Digestor

Product Code :
LBNY-0001-830037



Description :

COD Digestor

Technical Specification :

COD reactor is constructed of durable materials.
The aluminum block incorporates a 25-vial capacity and a well for a reference temperature probe.
The COD reactor is an easy to use test tube heater.
Its well-marked user interface provides intuitive operation.
The reactor is equipped with two predefined temperature settings.
Reactor Type; Dry temperature control device
Number of vial cells; 25 vial cell with 16mm holes
Programmable temp. range: RT +5°C -160°C
Temperature accuracy (°C): ±0.5°C at 40°, ±1°C at 120°C)
Temperature Stability; ±0.5°C
Heating time; ?15min (from 25°C to 160°C)
Display accuracy; ±1°C
Time range: 1 min ~ 99h59min
Heating Block; Aluminium
Weight (kg); 5.0 kg
Dimensions (mm); 260 x 220 x 100 (mm)
Power; 300W
Power Requirement; AC 220V, 50/60Hz.



Laboratory Instrument India