



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Digital Hot Plate Stirrer

**Product Code :**  
LBNY-0001-830038



**Description :**

Digital Hot Plate Stirrer

**Technical Specification :**

Digital is the perfect solution for temperature-sensitive laboratory applications where precise and accurate thermoregulation of the sample is required to ensure maximum reproducibility.

User interface: digital

Speed control: 1 digital

Stirring speed range: 50 - 1500 rpm

Stirring volume: 20 L

Temperature range: room temp to 310 °C

Temperature setting resolution: 1 °C

Thermoregulation accuracy:  $\pm 0.5$  °C (vtf);  $\pm 1.0$  °C (pt100, pt100)

Hot plate warning system: On Display.



**Laboratory Instrument India**

