



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Acid-resistant Gloves

**Product Code :**  
LBNY-0007-101000130



**Description :**

Rubber products, long sleeves with five fingers.  
The sleeve length is not less than 30cm corrosion resistant to acid, alkali and oxidant and reducing agent such as chemical reagent, and durable.

**Technical Specification :**

Acid-resistant Gloves



**Laboratory Instrument India**