



Email : sales@laboratoryinstrumentindia.com

Product Name :
Autoclave

Product Code :
LBNY-0001-830041



Description :

Autoclave

Technical Specification :

Made from heavy gauge stainless steel, this autoclave generates steam under high-pressure and temperature to ensure complete sterilization.

The insulated jacket provided with the autoclave minimizes temperature loss, making the process more effective and quick.

Steam sterilizer, Vertical loader

Lock lever

Caster: with stopper

With pressure gauge

Drain valve: manual value

Chamber capacity: 105-130L

Useable max. temperature & pressure: 137°C & 0.270 Mpa

Useable time: sterilization/Dissolution: 0min – 48hrs. 00min (& continuous), warming: 00min – 48hrs. 00min

Controller (OPERATION PANEL): LCD display Microprocessor control display temperature- 3- digit setting, 3 digit display, PID control, timer- 2 digits hours: 2 digits minutes setting and display, down count temperature/ time integration system

Safety Device: Lid lock system, Double interlock system, check system of lid close, over pressure detector, over temperature detector, water lack preventor, sensor break indicator, abnormal time indicator, memory of power failure, safety valve, breaker for slight leakage, over current & short circuit

Standard accessories: Stainless steel wire basket (2 pieces), flexible drain hose

Material of chamber: stainless steel

Material of enclosure: steel with melamin resin/ heat resistant plastic.



Laboratory Instrument India