



Email : sales@laboratoryinstrumentindia.com

Product Name :
Vertical Electrophoresis Tank

Product Code :
LBNY-0008-1010009



Description :

Polycarbonate injection molding.

Small glass board, 100*74mm 4 pieces.

Inclined plate 2, wire 1 pay.

Product size: 140 * 100 * 120mm, Gel plate specification (L * W): 105*185mm. Sample box: 15 teeth, 10 teeth, 4 handles each.

Technical Specification :

Vertical Electrophoresis Tank



Laboratory Instrument India