




Email : sales@laboratoryinstrumentindia.com

Product Name :
Horizontal Electrophoresis Tank

Product Code :
LBNY-0008-10100010



Description :

The design of the transparent top cover is opening type. Polycarbonate injection molding.
Gel plate specification (L * W): 60 * 60mm, respectively 6 teeth, 11 teeth each.
No leakage gel tray 2, wire 2, removable electrode rack and electrode head.
Product size: 190 * 100 * 60mm.

Technical Specification :

Horizontal Electrophoresis Tank



Laboratory Instrument India