




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Horizontal Electrophoresis Tank

**Product Code :**  
LBNY-0008-10100010



**Description :**

The design of the transparent top cover is opening type. Polycarbonate injection molding.  
Gel plate specification (L \* W): 60 \* 60mm, respectively 6 teeth, 11 teeth each.  
No leakage gel tray 2, wire 2, removable electrode rack and electrode head.  
Product size: 190 \* 100 \* 60mm.

**Technical Specification :**

Horizontal Electrophoresis Tank



**Laboratory Instrument India**