



Email : sales@laboratoryinstrumentindia.com

Product Name :
Light Incubator

Product Code :
LBNY-0008-10100015



Description :

It is used in school biological laboratory, plant germination, seedling raising, tissue culture, microorganism culture, rearing of insects and small animals, BOD determination of water body analysis, etc..
Light incubator is a high precision constant temperature equipment with illumination function.
Programmable preset temperature, light intensity, and can set 99 control periods, each set of time range of 1-99 hours, you can simulate the day and night temperature, light changes.
The product adopts microcomputer intelligent control.
Heating power: 1000W;
Light power: 48W;
The illumination ?2500Lx;
Temperature control range: full illumination 10?50? / no light 5?50?;
Temperature accuracy: $\pm 1?$;
Temperature uniformity: 2?;
Range of illumination: 0 - 12000Lx;
Compressor power: 75W.

Technical Specification :

Light Incubator

Laboratory Instrument India

