



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Light Incubator

**Product Code :**  
LBNY-0008-10100015



**Description :**

It is used in school biological laboratory, plant germination, seedling raising, tissue culture, microorganism culture, rearing of insects and small animals, BOD determination of water body analysis, etc..  
Light incubator is a high precision constant temperature equipment with illumination function.  
Programmable preset temperature, light intensity, and can set 99 control periods, each set of time range of 1-99 hours, you can simulate the day and night temperature, light changes.  
The product adopts microcomputer intelligent control.  
Heating power: 1000W;  
Light power: 48W;  
The illumination ?2500Lx;  
Temperature control range: full illumination 10?50? / no light 5?50?;  
Temperature accuracy:  $\pm 1$ ?;  
Temperature uniformity: 2?;  
Range of illumination: 0 - 12000Lx;  
Compressor power: 75W.

**Technical Specification :**

Light Incubator

---

# Laboratory Instrument India

