



Email : sales@laboratoryinstrumentindia.com

Product Name :
Field Microscope

Product Code :
LBNY-0001-830042



Description :

Field Microscope

Technical Specification :

Viewing Head: Binocular or Trinocular Siedentopf type, inclined 30°; interpupillary distance adjustment 48-75mm, 360° rotatable

Eyepieces: Plan 10x eyepiece, 20mm field of view with built-in diopter adjustment; a pointer is standard in one eyepiece; eyepieces are locked-in to prevent classroom tampering. Eyepiece accepts a Ø20.5mm reticle

Nosepiece: Rear facing quadruple, ball bearing nosepiece with positive click stops to ensure parcentration

Objectives 1X, 4X, 10x and 40X.

Cordless uses rechargeable LED illumination.

Easy to switch from compound to macro by changing objectives and stage position.

Head inclination: range 30 degrees to 45 degrees

LED light rated up to 50 000 hours of life.

All metal construct for rugged outdoor life.

Supplied with batteries and battery charger.



Laboratory Instrument India

