



Email : sales@laboratoryinstrumentindia.com

Product Name :
Inoculation Box

Product Code :
LBNY-0008-10100023



Description :

Size is not less than 400mm * 500mm * 350mm, wooden structure.
The front side of the box body operation hole two, and with a white cloth sleeve connection and sealing, the sleeve is not less than 300mm.
One-person operation case, open two round holes on the front and bring an ultraviolet lamp.
Operating voltage: AC220V 50Hz, UV sterilization tube is 15W.

Technical Specification :

Inoculation Box



Laboratory Instrument India