



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Inoculation Box

**Product Code :**  
LBNY-0008-10100023



**Description :**

Size is not less than 400mm \* 500mm \* 350mm, wooden structure.  
The front side of the box body operation hole two, and with a white cloth sleeve connection and sealing, the sleeve is not less than 300mm.  
One-person operation case, open two round holes on the front and bring an ultraviolet lamp.  
Operating voltage: AC220V 50Hz, UV sterilization tube is 15W.

**Technical Specification :**

Inoculation Box



**Laboratory Instrument India**