



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Automatic Kjeldahl Analyzer

**Product Code :**  
LBNY-0001-830043



**Description :**

Automatic Kjeldahl Analyzer

**Technical Specification :**

The fully automatic process ensures efficient operations, distillation and titration performed simultaneously  
Shortest time-to-results with online titration and automatic results calculation  
Premium result accuracy and precision thanks to the integrated colorimetric titrator with high precision burette  
Programmable boric acid, water, sodium hydroxide addition  
Kjeldahl nitrogen determination  
Measuring range 0.04-220 mg N  
Detection Limit: ? 0.015 mg N  
Display: color touch screen  
Analysis Time: from 4 minutes (including titration)  
Delay time: 00 seconds- 99 minutes 59 seconds  
Fully automatic reagents addition  
Steam flow regulation 10-100%  
Distillation and Titration RESIDUES REMOVAL: automatic.

---

# Laboratory Instrument India

