



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Infrared Thermometer With Laser

**Product Code :**  
LBNY-0002-910001



#### **Description :**

Infrared Thermometer With Laser, Infrared Thermometer With Laser Equipments, Infrared Thermometer With Laser Tools, Infrared Thermometer With Laser Tool Kits, Infrared Thermometer With Laser Manufacturers, Infrared Thermometer With Laser Suppliers from India, China, Kenya

Lab Tender Suppliers Africa and Lab Equipment Africa Kenya, Zambia, Chad, Somalia, Senegal, Zimbabwe, Guinea, Rwanda, Benin, Burundi, Tunisia, Uganda, Sudan, Algeria, Morocco, Angola, Ghana, Mozambique, Madagascar, South Africa, Ivory Coast, Cameroon, Niger, Burkina Faso, Mayotte, Sao Tome, Seychelles, South Sudan, Mali, Malawi, Nigeria, Ethiopia, Egypt, DR Congo, Mauritania, Eritrea, Gambia, Botswana, Namibia, Gabon, Lesotho, Guinea Bissau, Equatorial Guinea, Mauritius, Eswatini, Djibouti, Reunion, Comoros, Cape Verde, Western Sahara, Tanzania, Togo, Sierra Leone, Libya, Congo, Central African Republic, Liberia

#### **Technical Specification :**

Infrared thermometer with laser

Range: Large temperature measuring range switch between fahrenheit and celsius temperature scales.

Spot Ratio: The distance to spot ratio is 12:1, it can accurately measure targets at greater distances compared to most other ir thermometers.

Accuracy: Built-in laser for precisely aiming built in laser gives you the precision to hone in on the exact space you want to measure. operates on a twelve-to-one ratio ,accuracy +/- 2% on 95 emissive surfaces ,response time of 500ms.

Auto Shut Off: Big backlit screen makes it very easy to see the measured temperature, feature auto shut-off to save battery and data hold function.

Mode Button: Mode button to hold the max and min data. temperature readings from this devices are inanimate

---

objects. The measured temperature for humans or animals will not be correct.

Non-contact temperature range: -58 to 1022°F (-50 to 550°C)

Field of view: 12: 1

Display illumination: Back-lit LCD color display

Target indicator: Circular laser spot

Emissivity: 0.95

Resolution: 0.1°F / 0.1°C

Accuracy:  $\pm 4^\circ\text{F}$  ( $2^\circ\text{C}$ ) or  $\pm 2\%$  of reading

Response time: 0.5 seconds

Power off: Automatic after 15 seconds

Battery life: 14 hours

Included accessories: 2 AAA batteries, instruction manual, plastic molded case.



**Laboratory Instrument India**