



Email : sales@laboratoryinstrumentindia.com

Product Name :
Analog Multimeter

Product Code :
LBNY-0002-620016



Description :

Analog Multimeter

Technical Specification :

Measures AC, DC Voltage
Audible continuity test
Compact & Lightweight
Overload protection
Clear display with high contrast needle
Measurement ranges:
DC Voltage: 1000V
AC Voltage: 750V
Resistance: 1K?

Analog Multimeter, Analog Multimeter Equipments, Analog Multimeter Tools, Analog Multimeter Tool Kits, Analog Multimeter Manufacturers, Analog Multimeter Suppliers from India, China, Kenya



Laboratory Instrument India

