



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Temperature and Humidity Data Logger

**Product Code :**  
LBNY-0002-910003



**Description :**

Temperature and Humidity Data Logger

**Technical Specification :**

Temperature and humidity Data logger  
Status Indication via Red/Yellow LED and Green LED  
USB Interface for Setup and Data Download  
User Programmable Alarm Thresholds for Relative Humidity and Temperature  
Analysis software used to view graph for logging data  
Multimode to start logging  
Long battery life  
12 hour data logging duration  
Temperature range -30 to 50 °C  
Humidity range 0 to 100% RH  
Best with li-ion rechargeable Battery if not AAA or AA battery  
With Microsoft Exile data output option.  
Digital.

---

# Laboratory Instrument India

