



Email : sales@laboratoryinstrumentindia.com

Product Name :
Screen and Support

Product Code :
LBNY-0009-10100044



Description :

The instrument shelter cover is horizontal two wood inlaid together, above the front behind the 10mm, in order to protect the equipment from damage.

The size of the instrument shelter is about 460mm * 290mm * 537mm (height, width and depth).

The height of the frame is 1500mm, the width of the material is 35mm, the thickness of 3mm is high quality carbon structural steel, and the surface is rustproof.

Technical Specification :

Screen and Support



Laboratory Instrument India